



Container:	Image:	Example Applications:
1. Blood Culture		Blood Culture.
2. Serum Gel 4.7 mL		B12, Folate, Cardiac Enzymes, Hepatitis A, B, C, HIV, Lithium, Troponin, Amylase, BHCG, CHOL/TRIGS, Drug Levels, Fe Studies, Hormones, LFT, SUA, TFT, U+E's, CEA.
3. Sodium Citrate 3 mL		APTT, D Dimer, Fibrinogen, INR.
4. Lithium Heparin tube 7.5 mL		Troponin (Cobram and Echuca), Chromosomes.
5. Metals tube 6 mL		Lead, Copper, Aluminium, Zinc, Mercury. NOTE: Use metals needles available from GVH Pathology.
6. EDTA 2.6 mL		Direct Coombs Test, ESR, FBE, HBA1c, IM Test, Red Cell Folate, Homocysteine.
7. EDTA 7.5 mL		Antibody Screens, Blood Groups, Cross Matching.
8. Fluoride Oxalate 2.6 mL		Alcohol, Glucose Challenge Test, Glucose Tolerance Test.
9. ACD 9 mL		HLA Typing, Coeliac Genetic Test.

For further information please contact the GV Health Pathology Collection Centre on (03) 5832 2358.

Online Collection Manual: <http://pathology.gvhealth.org.au>